

NEW Product Announcement

Single Use Wide Bore Transfer Pipet

Designed for the rapid and effective transfer of viscous liquids and small particulates.



Globe Scientific is pleased to introduce these disposable wide bore transfer pipets. These pipets feature a wide tip opening designed specifically for the transfer of viscous liquids and small particulates. They are useful in the biofuel, petroleum, water treatment, food/beverage and pathology industries. The pipets are manufactured from low-density polyethylene (LDPE) and are packaged in convenient dispenser boxes; 500 pipets/box, 10 boxes/case.

Ideal for transferring:

- Algae
- Sludge, oil, slurry
- Small tissue biopsies and pleural specimens
- Plant specimens
- Sauces, gravies, syrup and condiments

Item #	Description	Packaging Options
135040	Wide bore transfer pipet, large bulb, 500/dispenser box, 10 boxes/case	500 & 5000

