

# NEW Product Announcement

## 15mL Large Capacity Transfer Pipet

*Specially designed for water testing and other procedures where a larger capacity liquid transfer is required.*



Globe Scientific is pleased to introduce these large capacity, graduated transfer pipets. These disposable pipets feature a **generous 5mL bulb draw** and are an economical alternative to expensive serological pipettes when extreme accuracy isn't critical. The pipets are made from low-density polyethylene (LDPE) and are packaged 250 pipets per box, 10 boxes per case.

- Ideal for water testing procedures
- Ideal alternative to more expensive serological pipets when extreme accuracy isn't critical
- 5mL bulb draw
- Graduated to 5mL
- Non-Sterile

Item #	Description	Packaging Options
139060	Transfer pipet, 15mL, large bulb, 250/box, 10 boxes/case	250 & 2500

### Technical Data

- Capacity: 15 mL
- Drops/mL: 21
- Bulb draw: 5 mL
- Total length: 216 mm (8.5")
- Stem diameter: 7.6 mm (0.30")
- Tip diameter: 3.1 mm (0.12")
- Graduated up to: 5 mL
- Graduation increments: 1 mL

Drawing is actual size



[www.globescientific.com](http://www.globescientific.com)

Globe Scientific Inc.

610 Winters Avenue • Paramus, NJ 07652

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail [mail@globescientific.com](mailto:mail@globescientific.com)