

# Product**CHANGE** Announcement

## Cardboard Storage Boxes

**Redesigned boxes now feature internal numbered grids to aid in sample identification and weep holes to facilitate drainage.**


**Our popular cardboard freezer boxes have been redesigned to incorporate these added features:**

- A numbered grid is now printed on the inside bottom of each box. This will help users identify sample vials by location inside the box.
- Weep holes help drain excess water condensation or liquid nitrogen, reduce moisture build-up and provide better air flow inside the box. This will help the boxes last longer.

These rugged cardboard boxes are designed for the compact storage of sample tubes, cryogenic vials and centrifuge tubes in a temperature range of -196° to +121°C. They have a liquid resistant coating and are designed for use in refrigerator and cryogenic freezer systems.

- Boxes come fully assembled
- Each box has a liquid resistant coating
- Temperature range: -196° to +121° C

Ordering information on next page . . .



*Numbered grid aids in sample identification.*



*Weep holes drain liquid and reduce moisture build-up resulting in longer box life.*

## Cardboard Storage Boxes

These rugged cardboard boxes are designed for the compact storage of sample tubes, cryogenic vials and centrifuge tubes in a temperature range of -196° to +121°C. They have a liquid resistant coating and are designed for use in refrigerator and cryogenic freezer systems.

- Boxes come fully assembled
- Each box has a liquid resistant coating
- Temperature range: -196° to +121° C
- Numbered grid printed in bottom of each box
- Weep holes for reduced moisture build-up

### Storage Box for 2" Tall Tubes

Dimensions: 134 x 134 x 47mm

Item #	Description	Units
3090	64 place (8x8), white, for up to 2" tall x 15mm wide tubes	Each & 96
3091	81 place (9x9), white, for up to 2" tall x 13mm wide tubes	Each & 96
3091R	81 place (9x9), red, for up to 2" tall x 13mm wide tubes	Each & 96
3092	100 place (10x10), white for up to 2" tall x 12mm wide tubes	Each & 96

### Storage Box for 3" Tall Tubes

Dimensions: 134 x 134 x 76mm

Item #	Description	Units
3094	64 place (8x8), white, for up to 3" tall x 15mm wide tubes	Each & 48
3094Y	64 place (8x8), yellow, for up to 3" tall x 15mm wide tubes	Each & 48
3095	81 place (9x9), white, for up to 3" tall x 13mm wide tubes	Each & 48
3096	100 place (10x10), white for up to 3" tall x 12mm wide tubes	Each & 48

### Storage Box for 15mL Centrifuge Tubes

Dimensions: 134 x 134 x 115mm

Item #	Description	Units
3098	49 place (7x7), white, for 15mL centrifuge tubes	Each & 36

### Storage Box for 50mL Centrifuge Tubes

Dimensions: 134 x 134 x 122mm

Item #	Description	Units
3099	16 place (4x4), white, for 50mL centrifuge tubes	Each & 48

Additional colors available by special order.



[www.globescientific.com](http://www.globescientific.com)

**Globe Scientific Inc.**  
610 Winters Avenue • Paramus, NJ 07652

**Customer service, sales and technical support:**

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail [mail@globescientific.com](mailto:mail@globescientific.com)