

Diamond RV-Pette-PRO™

120 dispensing volume options!

Dosage Chart

globe[®]
SCIENTIFIC INC.
www.globescientific.com

		PROGRAM NO.									
		1	2	3	4	5	6	7	8	9	10
TIP SIZE	ITEM #	VOLUME DISPENSED (µL)									
0.05 mL	3920 / 3920S	1	1.5	2	2.5	3	3.5	4	4.5	5	5.5
0.10 mL	3921 / 3921S	2	3	4	5	6	7	8	9	10	11
0.20 mL	3922 / 3922S	4	6	8	10	12	14	16	18	20	22
0.50 mL	3923 / 3923S	10	15	20	25	30	35	40	45	50	55
1.00 mL	3924 / 3924S	20	30	40	50	60	70	80	90	100	110
1.25 mL	3925 / 3925S	25	37.5	50	62.5	75	87.5	100	112.5	125	137.5
2.50 mL	3926 / 3926S	50	75	100	125	150	175	200	225	250	275
5.00 mL	3927 / 3927S	100	150	200	250	300	350	400	450	500	550
10.00 mL	3928 / 3928S	200	300	400	500	600	700	800	900	1000	1100
12.50 mL	3929 / 3929S	250	375	500	625	750	875	1000	1125	1250	1375
25.00 mL	3930 / 3930S	500	750	1000	1250	1500	1750	2000	2250	2500	2750
50.00 mL	3931 / 3931S	1000	1500	2000	2500	3000	3500	4000	4500	5000	5500
		48	35	23	19	15	13	11	9	8	7
		STEPS									