


Diamond RV-Pette-PRO™

Dispenser Syringe Tips

COMPATIBILITY CHART

Wide range of sizes match all sizes offered by other manufacturers and eliminates the need to change protocols or practice!



Diamond RV-Pette-PRO™ Dispenser Syringe Tip Size:

| Pipettor Brand/Model: | 0.05mL | 0.1mL | 0.2mL | 0.5mL | 1.0mL | 1.25mL | 2.5mL | 5.0mL | 10.0mL | 12.5mL | 25mL | 50mL |
|------------------------------------|--------|-------|-------|-------|-------|--------|-------|-------|--------|--------|------|------|
| Diamond RV-Pette™ | YES | YES | YES | YES | YES | YES | YES | YES | YES | YES | YES | YES |
| Diamond RV-Pette PRO™ | YES | YES | YES | YES | YES | YES | YES | YES | YES | YES | YES | YES |
| Bel-Art Roxy-M™ Repeating Pipettor | ✓ | — | — | ✓ | — | ✓ | — | ✓ | — | ✓ | ✓ | ✓ |
| Capp Rhythm™ Repeater | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Capp R10™ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Brand® Handystep® S | — | ✓ | — | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Brand® Handystep® Electronic | — | ✓* | — | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* |
| Biohit® Mechanical Stepper | ✓ | — | — | ✓ | — | ✓ | ✓ | ✓ | — | ✓ | ✓ | ✓ |
| eLine® Light Dispenser | — | ✓ | ✓ | ✓ | ✓ | — | ✓ | ✓ | ✓ | — | ✓ | ✓ |
| eLine® Pro Dispenser | — | ✓ | ✓ | ✓ | ✓ | — | ✓ | ✓ | ✓ | — | ✓ | ✓ |
| Eppendorf® Repeater 4780 | ✓ | — | — | ✓ | — | ✓ | ✓ | ✓ | — | ✓ | ✓ | ✓ |
| Fisherbrand™ | — | ✓ | — | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Gilson® Distriman® | — | — | — | — | — | ✓ | — | — | — | ✓ | — | — |
| Gilson® Repetman® | — | ✓* | — | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* |
| Nichiryo NichiMate | ✓ | — | — | ✓ | — | ✓ | ✓ | ✓ | — | ✓ | ✓ | ✓ |
| Rainin® AutoRep S | — | ✓ | — | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| Rainin® AutoRep E | — | ✓* | — | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* | ✓* |
| VWR® Repeating Dispenser | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |

*Tip-size needs to be entered manually.

All instruments and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only. Globe Scientific Inc. expressly disclaims any affiliation to them.

