

Technical Data Sheet


Rev. 16040

Page 1 of 1

16 x 92.5mm (12mL Capacity) Flat Bottom Polypropylene (PP) Tube

Item# 115060





Description

These Globe Scientific flat bottom, self-standing tubes are made from unbreakable polypropylene (PP). The tubes are autoclavable. They are non-graduated, without rim.

Features

- 12mL capacity
- 16 x 92.5mm
- Polypropylene (PP)
- Non-graduated
- No rim
- Self-standing

Technical Specifications

Autoclavable: Yes

Brittleness temperature: 0°C

Maximum temperature: 135°C

Sterility: Non-sterile

Products

Item#	Description	Unit
115060	12mL Flat bottom tube, PP	1000



Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406 E-mail mail@globescientific.com 400 Corporate Drive • Mahwah NJ 07430

www.globescientific.com