

Graduated Transfer Pipet, 3mL

Item# 137035

Technical Specifications

Material:	Low Density Polyethylene (LDPE)
Length:	143mm (+/- 0.5mm)
Capacity:	3mL
Drops/mL:	25
Bulb Draw:	2.2mL
Graduated Up To:	1mL
Graduation Increments:	0.25mL
Sterility:	Non-Sterile
Autoclavable:	No
Temperature Range:	-50°C to +75°C

Chemical Resistance

- Excellent resistance (no attack) to dilute and concentrated acids, alcohols, bases and esters
- Good resistance (minor attack) to aldehydes, ketones and vegetable oils
- Limited resistance (moderate attack suitable for short-term use only) to aliphatic and aromatic hydrocarbons, mineral oils, and oxidizing agents
- Poor resistance, and not recommended for use with halogenated hydrocarbons.

Quality Control

- This product is lot inspected and tested to ensure the highest level of quality.

Recommended Storage Conditions

- Store in cool dry place
- Avoid temperature variations and humidity

