

3375 Series

# Manual Pipet Controllers



3375 Series Manual Pipet Controllers

The Globe Scientific 3375 Series of Manual Pipet Controllers are available in your choice of five colors. They provide easy, accurate filling and dispensing of serological pipettes between 0.1 to 100mL. Squeeze the bulb to fill your pipet, then aspirate by moving the analog lever — slightly for slow, more for fast. Blow out is accomplished by pressing the button just below the bulb. Your 3375 is protected by an easily replaced internal 3µm hydrophobic filter.

**Features:**

- Ease of use
- Precise pipetting control
- Robust and light weight
- Compatible with most plastic and glass pipettes from 0.1 - 100mL
- 3µm Replaceable hydrophobic filter
- Easy cleaning and maintenance
- Available in 5 colors

**ORDERING INFORMATION**

Item #	Description	Unit
3375B	Manual Pipette Controller, Blue	Each
3375G	Manual Pipette Controller, Green	Each
3375LB	Manual Pipette Controller, Light Blue	Each
3375R	Manual Pipette Controller, Red	Each
3375Y	Manual Pipette Controller, Yellow	Each
3375-FIL	3.0µm Filter	Each