



diamond SureFlow™

Laboratory Bottle Top Dispensers

For the best in bottle top dispensers, choose Globe Scientific's Diamond® SureFlow™!



ALL NEW!

Bottle Top Dispensers

Diamond[®] SureFlow[™] Bottle Top Dispensers have been designed with several features to make the dispensing of your valuable laboratory reagents easier than ever. With your choice of standard or recirculating models, features include a unique piston design which allows cleaning and servicing without changing the calibration, and a springless valve which requires less maintenance, offers higher reliability and added chemical resistance. Dual graduations, volume lock and four additional bottle thread adapters are all standard. Recirculating models help eliminate the loss of your valuable reagents by priming without loss.

FEATURES INCLUDE:

- Fully autoclavable
- 2-Year warranty
- Locking feature to prevent spills
- Single tool recalibrations
- 360° Rotation to dispense in multiple vessels
- Ergonomic design for easy dispensing
- Simple operation
- Lower pressure for larger volumes helps reduce splashing
- Telescopic fill tube for bottles ranging in height from 170mm to 300mm




Recirculating valve helps eliminate reagent loss

Large dispenser tube window display






GBTD-R-100 –
10.0mL to 100.0mL
with recirculating valve



GBTD-05 – Standard 0.5mL to 5.0mL



Volume lock
ensures accurate
measurements



Each Diamond® SureFlow™ package includes:
Bottle top dispenser, 4 threaded bottle adapters,
telescopic fill tube & calibration tool

diamond **SureFlow**™

diamond SureFlow™

Bottle Top Dispensers



Bottle Top Dispensers

Item #	Volume Range	Increments	Unit
Standard			
GBTD-2.5	0.25mL to 2.5mL	0.05mL	Each
GBTD-05	0.5mL to 5.0mL	0.1mL	Each
GBTD-10	1.0mL to 10.0mL	0.2mL	Each
GBTD-25	2.5mL to 25.0mL	0.5mL	Each
GBTD-30	2.5mL to 30.0mL	0.5mL	Each
GBTD-50	5.0mL to 50.0mL	1.0mL	Each
GBTD-60	5.0mL to 60.0mL	1.0mL	Each
GBTD-100	10.0mL to 100.0mL	2.0mL	Each

With Recirculation Valve

GBTD-R-2.5	0.25mL to 2.5mL	0.05mL	Each
GBTD-R-05	0.5mL to 5.0mL	0.1mL	Each
GBTD-R-10	1.0mL to 10.0mL	0.2mL	Each
GBTD-R-25	2.5mL to 25.0mL	0.5mL	Each
GBTD-R-30	2.5mL to 30.0mL	0.5mL	Each
GBTD-R-50	5.0mL to 50.0mL	1.0mL	Each
GBTD-R-60	5.0mL to 60.0mL	1.0mL	Each
GBTD-R-100	10.0mL to 100.0mL	2.0mL	Each

Accessories

Item #	Description	Unit
GBTD-ADP-28	28mm Bottle Top Adaptor	Each
GBTD-ADP-38	38mm Bottle Top Adaptor	Each
GBTD-ADP-40	40mm Bottle Top Adaptor	Each
GBTD-ADP-45	45mm Bottle Top Adaptor	Each
GBTD-CALTOOL	Replacement Calibration Tool	Each
GBTD-TUBE-RECIRC	Replacement Recirculation Tube	Each
GBTD-TUBE-TELES	Replacement Telescopic Fill Tube, 2.5mL to 60mL	Each



www.globescientific.com

Globe Scientific Inc.

400 Corporate Drive • Mahwah, New Jersey 07430 • USA

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com