

GTR-3D

# Nutating Tube Rocker

## Quick Start Guide

### What's in the box:

- Tube Rocker
- Power Cable(s)
- Quick Start Guide

### OPERATION

- Simply switch the main power to the rocker ON or OFF. When the unit is turned ON, it will rock in a nutating motion at a fixed speed of 20rpm.



### IMPORTANT:

- The tube rocker must not be overloaded or damage may occur. Please do not load the unit in excess of 5.5lbs (2.5kg).

[www.globescientific.com](http://www.globescientific.com)

**globe**<sup>®</sup>  
SCIENTIFIC INC.

### Customer Service, Sales and Technical Support:

Phone: 1 (800) 394-4562 • 1 (201) 599-1400

Fax: 1 (201) 599-1406 • E-mail: [mail@globescientific.com](mailto:mail@globescientific.com)