

Labeling for Cryogenic and Cold Storage

Cold storage labels use a special adhesive that is designed to withstand extreme cold temperatures as low as -196°C (-311°F). Our cold storage labels can be printed or written on and will not bleed, smudge or smear. These labels are offered in a wide variety of sizes and colors on rolls or laser sheets-- perfect for distinguishing various specimens after prolonged periods of cold storage.


Features:

- Able to withstand temperatures from -196°C up to $+80^{\circ}\text{C}$
- Can be printed on with a laser printer
- Can be written on
- Available on rolls or sheets
- Full line of colors
- Various sizes offered
- Printing instructions and template included




Cryogenic and Cold Storage Labels


Label Sheets




LCS-D95W



LCS-D13RW



LCS-24X13W



LCS-33X13W

Dots for Tubes

Item #	Description	Quantity	Color
LCS-D95W	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	White
LCS-D95B	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Blue
LCS-D95G	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Green
LCS-D95N	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Orange
LCS-D95R	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Red
LCS-D95Y	9.5mm, For .5-1.5ml Tubes	20 Sheets, 192 Labels Per Sheet	Yellow
LCS-D95RW	9.5mm, For .5-1.5ml Tubes	20 Sheets, 960 Labels per Color	Assorted
LCS-D13W	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	White
LCS-D13B	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Blue
LCS-D13G	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Green
LCS-D13N	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Orange
LCS-D13R	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Red
LCS-D13Y	13mm, For 1.5-2ml Tubes	20 Sheets, 192 Labels Per Sheet	Yellow
LCS-D13RW	13mm, For 1.5-2ml Tubes	20 Sheets, 960 Labels Per Color	Assorted

Assorted colors contains Blue, Green, Violet, Red, and Yellow


Labels for Tubes

Item #	Description	Quantity	Color
LCS-24X13W	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	White
LCS-24X13B	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Blue
LCS-24X13G	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Green
LCS-24X13N	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Orange
LCS-24X13R	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Red
LCS-24X13Y	24x13mm, For .5 Tubes	20 Sheets, 119 Labels Per Sheet	Yellow
LCS-24X13RW	24x13mm, For .5 Tubes	20 Sheets, 595 Labels Per Color	Assorted
LCS-33X13W	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	White
LCS-33X13B	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Blue
LCS-33X13G	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Green
LCS-33X13N	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Orange
LCS-33X13R	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Red
LCS-33X13Y	33x13mm, For 1.5-2ml Tubes	20 Sheets, 85 Labels Per Sheet	Yellow
LCS-33X13RW	33x13mm, For 1.5-2ml Tubes	20 Sheets, 425 Labels per Color	Assorted

Assorted colors contains Blue, Green, Violet, Red, and Yellow

Cryogenic and Cold Storage Labels

Label Sheets




LCS-25X25RW

For Microplates

Item #	Description	Quantity	Color
LCS-38X6W	38x6mm	20 Sheets, 156 Labels Per Sheet	White
LCS-38X6B	38x6mm	20 Sheets, 156 Labels Per Sheet	Blue
LCS-38X6G	38x6mm	20 Sheets, 156 Labels Per Sheet	Green
LCS-38X6N	38x6mm	20 Sheets, 156 Labels Per Sheet	Orange
LCS-38X6R	38x6mm	20 Sheets, 156 Labels Per Sheet	Red
LCS-38X6Y	38x6mm	20 Sheets, 156 Labels Per Sheet	Yellow
LCS-38X6RW	38x6mm	20 Sheets, 780 Labels per Color	Assorted
LCS-67X5W	67x5mm	20 Sheets, 99 Labels Per Sheet	White


Assorted colors contains Blue, Green, Violet, Red, and Yellow



LCS-25X25W

For Microscope Slides

Item #	Description	Quantity	Color
LCS-25X25W	25x25mm	20 Sheets, 80 Labels Per Sheet	White



LCS-38X19W

For General Use


Item #	Description	Quantity	Color
LCS-38X19W	38x19mm	20 Sheets, 60 Labels Per Sheet	White
LCS-38X19B	38x19mm	20 Sheets, 60 Labels Per Sheet	Blue
LCS-38X19G	38x19mm	20 Sheets, 60 Labels Per Sheet	Green
LCS-38X19N	38x19mm	20 Sheets, 60 Labels Per Sheet	Orange
LCS-38X19R	38x19mm	20 Sheets, 60 Labels Per Sheet	Red
LCS-38X19Y	38x19mm	20 Sheets, 60 Labels Per Sheet	Yellow
LCS-38X19RW	38x19mm	20 Sheets, 300 Labels per Color	Assorted

Assorted colors contains Blue, Green, Violet, Red, and Yellow

Cryogenic and Cold Storage Labels

Label Sheets

For Cryogenic Vials


Item #	Description	Quantity	Color
LCS-43X19W	43x19mm	20 Sheets, 52 Labels Per Sheet,	White 
LCS-43X19B	43x19mm	20 Sheets, 52 Labels Per Sheet	Blue 
LCS-43X19G	43x19mm	20 Sheets, 52 Labels Per Sheet	Green 
LCS-43X19N	43x19mm	20 Sheets, 52 Labels Per Sheet	Orange 
LCS-43X19R	43x19mm	20 Sheets, 52 Labels Per Sheet	Red 
LCS-43X19Y	43x19mm	20 Sheets, 52 Labels Per Sheet	Yellow 
LCS-43X19RW	43x19mm	20 Sheets, 260 Labels Per Color	Assorted 

Assorted colors contains Blue, Green, Violet, Red, and Yellow


For Racks and Boxes

Item #	Description	Quantity	Color
LCS-67X25W	67x25mm	20 Sheets, 30 Labels Per Sheet	White 
LCS-67X25B	67x25mm	20 Sheets, 30 Labels Per Sheet	Blue 
LCS-67X25G	67x25mm	20 Sheets, 30 Labels Per Sheet	Green 
LCS-67X25N	67x25mm	20 Sheets, 30 Labels Per Sheet	Orange 
LCS-67X25R	67x25mm	20 Sheets, 30 Labels Per Sheet	Red 
LCS-67X25Y	67x25mm	20 Sheets, 30 Labels Per Sheet	Yellow 
LCS-67X25RW	67x25mm	20 Sheets, 150 Labels Per Color	Assorted 
LCS-28X44W	28x44mm	42 Labels Per Sheet	White 
LCS-64X13W	64x13mm	20 Sheets, 60 Labels Per Sheet	White 
LCS-76X51W	76x51mm	20 Sheets, 12 labels Per Sheet	White 

Assorted colors contains Blue, Green, Violet, Red, and Yellow



LCS-38X19RW



LCS-67X25RW

Cryogenic and Cold Storage Labels

Label Rolls



LCR-95W



LCR-95RW



LCR-13R



LCR-19W

LCR-95BIO



LCR-95RAD



Dots for Tubes

Item #	Description	Quantity	Color
LCR-95W	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	White
LCR-95B	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Blue
LCR-95C	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Coffee
LCR-95GL	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Gold
LCR-95G	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Green
LCR-95L	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Lavender
LCR-95N	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Orange
LCR-95P	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Pink
LCR-95R	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Red
LCR-95SL	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Silver
LCR-95Y	9.5mm, For .5-1.5ml Tubes	1000 Labels per Roll	Yellow
LCR-95RW	9.5mm, For .5-1.5ml Tubes	1000 Dots each color	Assorted
LCR-13W	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	White
LCR-13B	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Blue
LCR-13C	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Coffee
LCR-13GL	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Gold
LCR-13G	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Green
LCR-13L	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Lavender
LCR-13N	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Orange
LCR-13P	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Pink
LCR-13R	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Red
LCR-13SL	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Silver
LCR-13Y	13mm, For 1.5-2ml Tubes	1000 Labels per Roll	Yellow
LCR-13RW	13mm, For 1.5-2ml Tubes	1000 Dots each color	Assorted
LCR-19W	19mm Dots, For 15-50ml Tubes	1000 Labels per Roll	White
LCR-21X7W	21x7mm, For 0.5ml Tubes	1500 Labels per Roll	White
LCR-24X13W	24x13mm, For 0.5ml Tubes	1000 Labels per Roll	White
LCR-25W	25mm Dots, For 50ml Tubes	500 Labels per Roll	White

Assorted colors contains Blue, Green, Violet, Red, and Yellow







Hazard Dots

Item #	Description	Quantity	Color
LCR-95BIO	9.5mm, Biohazard Dots	1000 Labels per Roll	Red/Black
LCR-95RAD	9.5mm, Radioactive Dots	1000 Labels per Roll	Yellow/Pink

Cryogenic and Cold Storage Labels

Label Rolls

For Tubes

Item #	Description	Quantity	Color
LCR-33X13W	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	White 
LCR-33X13B	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Blue 
LCR-33X13G	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Green 
LCR-33X13N	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Orange 
LCR-33X13R	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Red 
LCR-33X13Y	33x13mm, For 1.5/2ml Tubes	1000 Labels per Roll	Yellow 

LCR-33X13W



For General Use

Item #	Description	Quantity	Color
LCR-38X19W	38x19mm	1000 Labels per Roll	White 
LCR-38X19B	38x19mm	1000 Labels per Roll	Blue 
LCR-38X19G	38x19mm	1000 Labels per Roll	Green 
LCR-38X19N	38x19mm	1000 Labels per Roll	Orange 
LCR-38X19R	38x19mm	1000 Labels per Roll	Red 
LCR-38X19Y	38x19mm	1000 Labels per Roll	Yellow 

LCR-38X19R



Label Dispenser

Item #	Description	Quantity	Color
LCR-DIS	Multi-Roll Dispenser	1 each	Clear

LCR-DIS



Perfectly fits boxes of label rolls for easy dispensing and organization.



www.globescientific.com

Customer Service, Sales and Technical Support:

Phone: 1(800) 394-4562 · 1(201) 599-1400

Fax: 1(201) 599-1406


Email: mail@globescientific.com

Cryogenic and Cold Storage Labels

Direct Thermal Label Rolls

Direct thermal Cryogenic Storage Labels share the same features and benefits as our standard cryo line with one major added perk--they work in your Direct Thermal Printer!


With a wide variety of sizes, colors and configurations you are sure to find a solution perfect for your existing setup.




LTC-95W



LTC-38X13W



LTC-51X6W



LTC-22X22W

Dots for Tubes

Item #	Description	Quantity	Color
LTC-95W	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	White
LTC-95B	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Blue
LTC-95G	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Green
LTC-95N	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Orange
LTC-95R	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Red
LTC-95Y	9.5mm, For .5-1.5ml Tubes	4000 Labels per Roll	Yellow
LTC-13W	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	White
LTC-13B	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Blue
LTC-13G	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Green
LTC-13N	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Orange
LTC-13R	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Red
LTC-13Y	13mm, For 1.5-2ml Tubes	4000 Labels per Roll	Yellow

Labels for Cryogenic Vials

Item #	Description	Quantity	Color
LTC-27X13W	27x13mm	1000 Labels per Roll	White
LTC-27X13B	27x13mm	1000 Labels per Roll	Blue
LTC-27X13G	27x13mm	1000 Labels per Roll	Green
LTC-27X13N	27x13mm	1000 Labels per Roll	Orange
LTC-27X13R	27x13mm	1000 Labels per Roll	Red
LTC-27X13Y	27x13mm	1000 Labels per Roll	Yellow
LTC-38X13W	38x13mm	750 Labels per Roll	White
LTC-38X13B	38x13mm	750 Labels per Roll	Blue
LTC-38X13G	38x13mm	750 Labels per Roll	Green
LTC-38X13N	38x13mm	750 Labels per Roll	Orange
LTC-38X13R	38x13mm	750 Labels per Roll	Red
LTC-38X13Y	38x13mm	750 Labels per Roll	Yellow

Labels for Microplates

Item #	Description	Quantity	Color
LTC-51X6W	51x6mm	1000 Labels per Roll	White

Labels for Microscope Slides

Item #	Description	Quantity	Color
LTC-22X22W	22x22mm	1000 Labels per Roll	White














Cryogenic and Cold Storage Labels

Direct Thermal Label Rolls












LTC-38X19W

Labels for Large Tubes, Racks & Boxes

Item #	Description	Quantity	Color
LTC-38X19W	38x19mm	750 Labels per Roll	White 
LTC-38X19B	38x19mm	750 Labels per Roll	Blue 
LTC-38X19G	38x19mm	750 Labels per Roll	Green 
LTC-38X19N	38x19mm	750 Labels per Roll	Orange 
LTC-38X19R	38x19mm	750 Labels per Roll	Red 
LTC-38X19Y	38x19mm	750 Labels per Roll	Yellow 
LTC-43X19W	43x19mm	500 Labels per Roll	White 
LTC-51X19W	51x19mm	500 Labels per Roll	White 
LTC-51X19B	51x19mm	500 Labels per Roll	Blue 
LTC-51X19G	51x19mm	500 Labels per Roll	Green 
LTC-51X19N	51x19mm	500 Labels per Roll	Orange 
LTC-51X19R	51x19mm	500 Labels per Roll	Red 
LTC-51X19Y	51x19mm	500 Labels per Roll	Yellow 

Label/Dot Combo for Tubes

Item #	Description	Quantity	Color
LTC-24X13-95W	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	White  
LTC-24X13-95B	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Blue  
LTC-24X13-95G	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Green  
LTC-24X13-95N	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Orange  
LTC-24X13-95R	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Red  
LTC-24X13-95Y	24x13mm & 9.5mm Dot, For .5-1.5ml Tubes	500 Labels per Roll	Yellow  
LTC-38X13-95W	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	White  
LTC-38X13-95B	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Blue  
LTC-38X13-95G	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Green  
LTC-38X13-95N	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Orange  
LTC-38X13-95R	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Red  
LTC-38X13-95Y	38x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Yellow  
LTC-33X13-95W	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	White  
LTC-33X13-95B	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Blue  
LTC-33X13-95G	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Green  
LTC-33X13-95N	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Orange  
LTC-33X13-95R	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Red  
LTC-33X13-95Y	33x13mm & 9.5mm Dot, For 1.5-2ml Tubes	500 Labels per Roll	Yellow  



LTC-33X13-95W



www.globescientific.com

Customer Service, Sales and Technical Support:

Phone: 1(800) 394-4562 · 1(201) 599-1400

Fax: 1(201) 599-1406

Email: mail@globescientific.com