

NEW PRODUCT

NEW!

DiamondPerfectMark™

The Ideal Laboratory Marker for Fine Labeling on Small Surfaces.




The Diamond PerfectMark™ is a permanent laboratory marker specially suited for laboratory marking applications. The fine tip makes it easy to write on very small surfaces like histology cassettes, microscope slides and test tubes. The marker is smear-proof, waterproof, odorless and light-resistant. The permanent ink is resistant to most laboratory chemicals & solvents.

Features:

- Extra fine 0.75mm nib is perfect for clear labeling on small surfaces.
- Works well on smooth and transparent surfaces like microscope slides, test tubes and labware
- Quick-drying, smudge-proof, waterproof, lightfast and low-odor.
- Resistant to abrasion and high temperatures
- Black ink

Perfect for Histology & Cytology Applications



Item #	Description	Unit
1000	Diamond PerfectMark Laboratory Marker, Black	10/Box



globe®
SCIENTIFIC INC.

www.globescientific.com

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com