

NEW PRODUCT



NEW!

TransferTop™

Urine Collection and Direct Transfer Cup

TransferTop™ Urine Collection and Direct Transfer Cups feature an integrated transfer mechanism for “closed system” transfer of urine samples directly from the cup to a vacuum sample tube. With the TransferTop™ System, there is no need for laboratory personnel to open the collection cup and pour the sample into a tube for processing. Instead, the safety label is peeled back, and the vacuum tube inserted for simple, safe, efficient, and reliable sample transfer.

TransferTop™ Closed System for Urine Transfer:

- Requires no open transferring of samples for reduced spills, improved lab hygiene, and limited exposure to dangerous samples
- Reduces risk of contamination for improved test results
- Directly transfers samples to a tube suitable for analysis and transport
- Works with most vacuum urine tubes
- Provides consistent sample transfer volumes when used with like vacuum tubes
- Improves efficiencies and sample throughput when compared to traditional open systems

Item #	Description	Unit
3852	Urine Collection & Direct Transfer Cup, 2oz (60mL) capacity, STERILE, individually wrapped	500/Case
3856L	Urine Collection & Direct Transfer Cup, 4oz (120mL) capacity, STERILE, bulk	300/Case
3841	Urine Transfer Straw, 3" (7.5cm) straw, non-sterile, 100/bag, 8 bags/case	800/Case



3856L – 4oz (120mL) TransferTop™ Collection and Transfer Cup with integrated transfer mechanism (graduated to 100mL)



3852 – 2oz (60mL) TransferTop™ Collection and Transfer Cup



3841 – Urine Transfer Straw for closed transfer of urine when using a standard collection container



www.globescientific.com

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com