


SHORTY SEROLOGICAL PIPETTES



NEW!

GET SHORTY!

Shorty serological pipettes for improved ergonomics and enhanced freedom of movement when used in fume hoods, biosafety cabinets, isolators and other confined spaces.



Globe Scientific shorty serological pipettes are approximately 1/3 shorter than standard serological pipettes of equal volume!

GET SHORTY!

*Standard length serological pipettes
345mm long*

Shorty serological pipettes — 230mm long



1880 – 5mL

1884 – 10mL

1888 – 25mL

Shorty serological pipettes for improved ergonomics and enhanced freedom of movement when used in confined spaces.

Globe Scientific shorty serological pipettes are approximately 1/3 shorter than standard serological pipettes of equal volume. The pipettes are manufactured under strict ISO 9001:2008 standards in a class 100,000 cleanroom. They are non-pyrogenic, non-cytotoxic, RNase-free and DNase-free. The pipettes are equipped with an aerosol barrier plug that prevents both liquid and aerosols from contaminating the pipettor. Printed black graduations prevent ambiguity and guarantee precise reading with less than 2% inaccuracy at full volume. Calibration is certified for accuracy. Available in 5mL, 10mL and 25mL sizes.

Features at a Glance:

- All 230mm long
- Material: USP Class VI approved polystyrene (PS)
- Manufactured under strict ISO 9001:2008 standards
- Produced in a class 100,000 cleanroom
- Non-pyrogenic and non-cytotoxic
- RNase-free and DNase-free
- Sterile pipettes are sterilized by gamma irradiation (SAL 10⁻⁶)
- Aerosol barrier filter plug: hydrophilic PET
- Reverse graduations on all 3 sizes
- Negative graduations on all 3 sizes
- Calibrated to deliver (TD) — calibration is certified for accuracy
- Accuracy: ±2% at full volume
- Pipettes are designed with anti-drip tips
- Fits all major brands of pipette controllers
- Color coded



Packaging:

- Available sterile in polyethylene and medical paper single peel packs.
- Designed for easy peel apart or pop-through pipette removal.

Item #	Description	Unit
1880	5mL, PS, sterile, individually wrapped	250
1884	10mL, PS, sterile, individually wrapped	200
1888	25mL, PS, sterile, individually wrapped	100



Globe Scientific Inc.
610 Winters Avenue
Paramus, New Jersey 07652 • USA

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400
Fax 1 (201) 599-1406
E-mail mail@globescientific.com

www.globescientific.com