

# Translucent Snap Caps

for evacuated blood drawing tubes, glass and plastic tubes.

These colorful snap-on closures are an effective way to recap blood drawing tubes, glass test tubes and plastic culture tubes. They are precision molded from low-density polyethylene (LDPE) and maintain a tight seal during centrifugation, agitation and refrigeration procedures. They are available in bright translucent colors for sample identification purposes. The caps feature a single thumb tab for easy removal and close nesting in tube racks.

## Features:

- Economical way to recap glass tubes, regardless of manufacturer
- Prevents spills and evaporation by employing a liquid tight seal
- Available in seven translucent colors that correspond to Vacutainer<sup>®</sup> cap colors: Natural, Blue, Green, Lavender, Red, Gray, and Yellow
- Caps can be written on with indelible ink for easy sample identification



Single thumb tab allows the tubes to be nested closely in test tube racks and still provides easy removal of the cap.



**See next page for Product Selection Guide >>**

**Translucent Snap Caps**  
**Product Selection Guide**



These caps			Fit these tube types and sizes		
Item #	Color	Cap Size	Evacuated Tube Size	Glass Culture Tube Size	Plastic Tube Size
113142A	Gray	12mm	13mm	12mm	12mm
113142B	Blue	12mm	13mm	12mm	12mm
113142C	Natural	12mm	13mm	12mm	12mm
113142G	Green	12mm	13mm	12mm	12mm
113142L	Lavender	12mm	13mm	12mm	12mm
113142R	Red	12mm	13mm	12mm	12mm
113142Y	Yellow	12mm	13mm	12mm	12mm
<hr/>					
113143A	Gray	13mm	—	13mm	—
113143B	Blue	13mm	—	13mm	—
113143C	Natural	13mm	—	13mm	—
113143G	Green	13mm	—	13mm	—
113143L	Lavender	13mm	—	13mm	—
113143R	Red	13mm	—	13mm	—
113143Y	Yellow	13mm	—	13mm	—
<hr/>					
113144A	Gray	16mm	16mm	—	—
113144B	Blue	16mm	16mm	—	—
113144C	Natural	16mm	16mm	—	—
113144G	Green	16mm	16mm	—	—
113144L	Lavender	16mm	16mm	—	—
113144R	Red	16mm	16mm	—	—
113144Y	Yellow	16mm	16mm	—	—



All caps sold in units of 1000 caps per bag.