


16 x 100mm (12mL) Polystyrene (PS) Test Tube

Item# 111012



Description

These Globe Scientific test tubes are molded from premium grade polystyrene (PS) for clarity and can withstand 1500 x g during centrifugation.

Features

- 12mL capacity
- 16 x 100mm
- Non-graduated
- No rim

Technical Specifications

Material: Polystyrene (PS)
 Spin Rate: 1500 x g
 Autoclavable: No
 Brittleness temperature: 20°C
 Maximum Temperature: 90°C
 Sterility: Non-sterile

Products

| Item# | Description | Unit |
|--------|--------------------------------|------|
| 111012 | 16 x 100mm, 12mL test tube, PS | 2000 |