

16 x 100mm (12mL Capacity) Polypropylene Test Tube

Item# 111041



Description


These Globe Scientific test tubes are molded from premium grade polypropylene (PP) for impact strength and stability to common acids, alkalies, and solvents. They can withstand 3000 RPM and are autoclavable and freezable.

Features

- 12mL capacity
- 16 x 100mm
- Polypropylene (PP)
- Non-graduated, no rim

Technical Specifications

Spin Rate: 3000 RPM
 Autoclavable: Yes
 Brittleness temperature: 0°C
 Maximum temperature: 135°C
 Sterility: Non-sterile



Products

| Item# | Description | Unit |
|--------|--------------------|------|
| 111041 | 12mL test tube, PP | 2000 |