

200-1000µL Pipette Tips for MLA[®], Natural

Item# 151147



151147

Description

These 200-1000 µL pipette tips are designed for use with MLA[®] pipettors. They are bulk packaged in bags of 1000 tips

Technical Specifications

Material: Virgin Polypropylene (PP)
 Autoclavable: Yes
 Temperature Range: -10°C - 120°C
 Graduations: Non-graduated
 Color: Natural
 Maximum Volume: 1000µL
 Sterility: Non-Sterile

Products

Item #	Description	Unit
151147	Pipette tip for MLA [®] pipettors, natural, 1000/pack	1000

