

## 2mL Nesting Sample Cup, Polystyrene (PS), for 16mm Tubes

Item# 5505



### Description

This nesting sample cup is used in conjunction with other test tubes or original blood collection tubes when smaller sampling is necessary.

It is ideal for use on direct-tube sampling machines: Olympus<sup>®</sup> AU480 and AU680, Dade Dimension<sup>®</sup>, Hitachi<sup>®</sup> 747, Kodak Ektachem<sup>®</sup> 700 Series and more.

### Features

- 2mL capacity
- For use with 16mm tubes
- Excellent clarity

### Products

Item#	Description	Unit
5505	2mL nesting sample cup, PS	1000

### Technical Specifications

Material: Polystyrene (PS)  
 Autoclavable: No  
 Temperature range: -10°C to +50°C  
 Sterility: Non-sterile

