

Reaction Tube for Sysmex[®] CA[®] Series Analyzers

Item# 5530



Description

This Globe Scientific polystyrene (PS) reaction tube is designed for use on the Sysmex[®] CA50, CA500, CA1000, CA1500 and CA5000 analyzers.

Features

- 0.8mL Capacity
- Excellent clarity

Technical Specifications

Material: Polystyrene (PS)
Dimensions: 10mm x 30mm
Capacity: 0.8mL
Autoclavable: No
Sterility: Non-sterile

Products

Item#	Description	Unit
5530	Reaction Tube, Polystyrene (PS), 1000 per bag	1000

Note: All instruments and names of manufacturers mentioned herein are the trademarks of their respective companies and are used for reference purposes only. Globe Scientific Inc. expressly disclaims any affiliation to them.