

## 0.5mL Sample Cup, Polystyrene (PS)

Item# 5541



### Description

These Globe Scientific multi-purpose sample cups are designed for use on the most popular instruments in the industry. The sample cups are made to exact specifications and provide excellent optical properties.

### Features

- 0.5mL capacity
- Excellent clarity

### Technical Specifications

Material: Polystyrene (PS)

Dimensions: 13.55mm x 24.50mm

Autoclavable: No

Brittleness temperature: 20°C

Maximum temperature: 90°C

Sterility: Non-sterile

### Products

Item#	Description	Unit
5541	0.5mL sample cup, PS	1000

