

## Fine Tip Transfer Pipet, 1.3mL

Item# 134010

### Technical Specifications

Material: Low Density Polyethylene (LDPE)  
 Length: 51mm (2.1")  
 Capacity: 1.3mL  
 Drops/mL: 45  
 Bulb Draw: 1mL  
 Sterility: Non-Sterile  
 Autoclavable: No  
 Temperature Range: -50°C to +75°C  
 Graduations: Non-Graduated

### Chemical Resistance


- Excellent resistance (no attack) to dilute and concentrated acids, alcohols, bases and esters
- Good resistance (minor attack) to aldehydes, ketones and vegetable oils
- Limited resistance (moderate attack suitable for short-term use only) to aliphatic and aromatic hydrocarbons, mineral oils, and oxidizing agents
- Poor resistance, and not recommended for use with halogenated hydrocarbons.

### Quality Control

- This product is lot inspected and tested to ensure the highest level of quality.

### Recommended Storage Conditions

- Store in cool dry place
- Avoid temperature variations and humidity



### Products

Item#	Description	Units
134010	Transfer pipet, 1.3mL, fine tip, 51mm, 1mL bulb draw	500 & 4000