

## Nexius Zoom Microscope Series

Item# ENZ-1703-S

### Eyepieces

- Supplied with a pair of HWF10x/23mm eyepieces.

### Head

- Trinocular head with 45° inclined tubes. Both eyepieces with +/-5 diopter adjustments. Interpupillary distance adjustable between 54mm and 75mm. Fixed light patch beamsplitter (50:50). The photo port is of Ø23.2mm

### Objectives

- Supplied with 1:8.4 zoom objective with 0.65x to 5.5x providing a total magnification of 6,5x until 55x, with a field of view from 35.4mm to 4.2mm and working distance 110mm.

*All optics are anti-fungus treated and anti-reflection coated for maximum light throughput.*

### Stand

- Ergonomic designed rack & pinion stand with 3W incident and 3W transmitted LED illuminations (-S stand)

### Illumination

- The 3W LED illuminators operate on an internal power supply of 100-240VAC. Both illuminations can be used simultaneously and the light intensities can be adjusted separately.

### Package Content

- Supplied with power cord, dust cover, a spare fuse and user manual. All packed in a polystyrene box.

