

16 x 101mm Transport Tube with Red Screw Cap, Sterile, Round Bottom, Self-Standing

Item# 6102S-R

Description

Sterile 16 x 101mm polypropylene (PP) transport tube with separate high-density polyethylene (HDPE) red screw cap. These leak-resistant tubes are tested to 70 Kpa and gamma sterilized. The self-standing tubes feature rounded bottoms with skirt.

Features

- Working volume 10mL to 13mL
- Brim Capacity 15mL
- Leak resistant – tested to 70 Kpa
- Sterilized by Irradiation
- Self-standing
- Non-graduated

Material Specifications:

- Tubes: USP VI Class approved polypropylene (PP)
Caps: High-density polyethylene (HDPE)
- DNase and RNase free
- Non-Pyrogenic

Quality Control:

- ISO 9001:2015 , ISO 13485:20 16

Packaging:

Tubes: 200 per heat-sealed individual bag, 5 bags per case

Caps: 1000 per bag, 1 bag per case

Products

Item #	Description	Unit
6102S-R	16 x 101mm transport tube, separate red screwcap, STERILE, self-standing, 200 tubes/bag, 1000 caps/case	1000

