

**NEW!**

# Four-Channel Laboratory Alarm Timer

The ideal timing solution for busy labs.



**Run up to four tests simultaneously!** This convenient timer permits setting 4 separate channels simultaneously in any combination of count down or count up (stopwatch). Set any count down time from 99 hours, 59 minutes, 59 seconds to 1 second with an accuracy of 0.01%. When zero is reached, the alarm sounds for one minute and timer starts counting up to show time elapsed since alarming. Each of the 4 channels has a distinctive alarm that can be turned off with a simple touch of any button. Additional functions include stopwatch, time-out, and time-of-day clock. Large keys make it easy to set times or change channels, and the 3/4" LCD display facilitates viewing from long distances.

## Features:

- Set up to 4 channels simultaneously
- Memory returns display to previously programmed count down time with the touch of a button
- Loud, 70 db alarm sounds for 1 minute or can be silenced manually
- Large LCD display allows viewing from across the lab
- 6 digit display shows hour, minute and second setting for timers and clock
- Choose between 12 or 24 hour clock format
- Provided with fixed stand, pocket clip, hole for lanyard, and magnetic back
- 1.5V G-13 Battery included



## Specifications

Channels/Zones	4
Timing Accuracy	±0.5 sec/day
Timing Resolution	1 second
Maximum Setting	99 hours, 59 minutes, 59 seconds
Clock Mode	12 and 24 hour format
Display	LCD
Material	ABS Plastic
Battery Type	1.5V G-13 size button cell battery (included)
Dimensions	2.75 in x 2.625 in x 1 in

## Functions

Memory
Count Up/Down
Time In/Out
Alarm
Time of Day (12 and 24 Hour Clock)

Item #	Description	Unit
1004	4 Channel Count Up-Count Down Timer with pocket clip, stand, magnetic holder, and battery	Each