

Sterile Low Retention

Filter Pipette Tips

All NEW!

150900

150903

150910

150915

150920

150925

150935

150938

These Globe Scientific STERILE filter pipette tips are manufactured from clear high-quality polypropylene (PP). Special processing renders these tips hydrophobic for exceptionally low sample retention. The tips are certified DNase free, RNase free and non-pyrogenic and are manufactured in a Class 100,000 cleanroom environment. A universal fit design provides an airtight fit on most brands of pipettors.

Features:

- Manufactured from clear polypropylene (PP)
- Polypropylene (PP) filters protect samples and pipettors from contamination
- DNase-free, RNase-free and non-pyrogenic
- Manufactured in a Class 100,000 cleanroom environment under ISO13485:2016 and ISO9001:2015 quality management systems
- Sterilized by gamma irradiation SAL10⁻⁶ (ISO11137)
- Racks are printed with lot number for quality traceability

The tips are offered in extra-rigid autoclavable racks and packaged 96/Rack, 10 Racks/Box, 2 Boxes/Carton



Item #	Description	Unit
150900	Filter Pipette Tip, 0.1-10uL, Certified, Universal, Low Retention, 32mm, Natural, STERILE	1920
150903	Filter Pipette Tip, 0.1-10uL Long, Certified, Universal, Low Retention, 45mm, Extended Length, STERILE	1920
150910	Filter Pipette Tip, 1-20uL, Certified, Universal, Low Retention, 50mm, Natural, STERILE	1920
150915	Filter Pipette Tip, 1-100uL, Certified, Universal, Low Retention, 50mm, Natural, STERILE	1920
150920	Filter Pipette Tip, 1-200uL, Certified, Universal, Low Retention, 59mm, Natural, STERILE	1920
150925	Filter Pipette Tip, 1-300uL, Certified, Universal, Low Retention, 59mm, Natural, STERILE	1920
150935	Filter Pipette Tip, 1-1000uL, Certified, Universal, Low Retention, 86mm, Natural, STERILE	1920
150938	Filter Pipette Tip, 1-1000uL Long, Certified, Universal, Low Retention, 105mm, Extended Length, Natural, STERILE	1920