

# Product **Spec Sheet**

**Globe Item#:** V-6297T

**Description:**


Centrifuge Tube, 50mL, Attached Blue Flat Top Screw Cap, Printed Graduations

**Material:** Polypropylene (PP), USP Class VI

**Specifications:**

Gamma Sterile  
RNase/DNase Free  
Cytotoxin Free, Non-Pyrogenic  
Temperature Range: -40°C to 121°C  
Max RCF: 16,000 x g

**Packaging:** 500/Case  
25/Bag, 20 Bags/Case



**Customer service, sales and technical support:**

Phone 1 (800) 394-4562 • 1 (201) 599-1400  
Fax 1 (201) 599-1406  
E-mail [mail@globescientific.com](mailto:mail@globescientific.com)

**Globe Scientific Inc.**

400 Corporate Drive • Mahwah, NJ 07430  
[www.globescientific.com](http://www.globescientific.com)